

PARTNER University of California, San Francisco (UCSF), Sophie Huddart, PhD, MSc	THERAPEUTIC AREA Tuberculosis / respiratory infectious disease
STUDY TYPE Multi-country, decentralized, prospective continuous cough monitoring, nested in a TB diagnostics trial (R2D2 TB Network)	PUBLICATION Peer-reviewed, published in the International Journal of Tuberculosis and Lung Disease
SETTING Five countries — Uganda, South Africa, Vietnam, the Philippines, India	YEAR Monitoring completed 2022; published 2023 (Int J Tuberc Lung Dis)

KEY OUTCOME

Across 565 participants in five countries, objective cough monitoring showed that TB cough falls within the first days of treatment, distinct from other respiratory disease, establishing continuous cough frequency as a fast, decentralized readout of TB treatment response.

5

COUNTRIES

565

PARTICIPANTS
MONITORED

23.9h

MEDIAN HOURS
MONITORED

14 days

MONITORING PERIOD

CHALLENGE

Cough is the entry point to almost every TB workup, yet it is one of the least reliable signals in respiratory medicine. Patients routinely misjudge how much they cough, and a substantial share of people with confirmed pulmonary TB report no cough at all. The tools used to confirm treatment response are objective but slow: sputum culture, the reference standard, can take six to eight weeks to return a result.

For a research team working across high-burden, low-resource settings, that left a gap. There was no lean, objective, decentralized way to measure whether and how fast a patient's cough was actually changing after treatment began, and no way to compare that trajectory between people who turned out to have TB and people who did not.

The UCSF team set out to test whether passive, continuous cough monitoring could fill that gap, in the real-world conditions where TB is actually diagnosed and treated rather than in a controlled clinic.

SOLUTION

Hyfe deployed its Research application on study-issued smartphones across five countries, nested inside the R2D2 TB Network's ongoing diagnostics trial. Adults presenting as TB suspects, with new or worsening cough of at least two weeks, were enrolled and monitored continuously for 14 days from the point of enrollment, spanning the initiation of TB treatment.


The application captured coughs passively as participants went about daily life, using an on-device two-stage algorithm that detects explosive sounds and classifies them as cough. Only cough timestamps left the device; no audio was transmitted, preserving participant privacy by design. TB status was established independently against a microbiological reference standard (sputum Xpert MTB/RIF Ultra and culture), allowing the team to compare cough trajectories between participants with microbiologically confirmed TB, clinically diagnosed TB, and other respiratory disease.

What made the deployment distinctive was that it worked under these conditions: a fully decentralized, multi-country, continuous-monitoring setup running in exactly the settings where TB care happens, generating objective, timestamped cough data at population scale.

RESULTS

Metric	Result
Participants monitored (with reference-standard result)	565
Countries	5
Continuous monitoring period	14 days
Microbiologically confirmed TB	144 (25.5%)
Clinically diagnosed TB	48 (8.5%)
Other respiratory disease	373 (66.0%)
Day 1 median cough/hour – confirmed TB	8
Day 1 median cough/hour – other respiratory disease	5
Day 1 difference (confirmed TB vs other respiratory disease)	$p < 0.001$
Overall median cough/hour, Day 1 → Day 14	5.0 → 3.5
Decline by Day 14 (all groups)	Significant ($p < 0.001$)

Headline finding: Participants with microbiologically confirmed TB began with markedly higher cough frequency and improved more rapidly once treatment started, following a trajectory distinct from other respiratory disease and converging with it by Day 14.



OPERATIONAL SNAPSHOT

MONITORING

14 days continuous, spanning treatment initiation

DESIGN

Decentralized, prospective, continuous cough monitoring nested in a multi-country TB diagnostics trial (R2D2 TB Network)

PARTICIPANTS

565 adults presenting as TB suspects (cough $\geq$ 2 weeks)

SITES

Five countries — Uganda, South Africa, Vietnam, the Philippines, India

TECHNOLOGY

Hyfe Research application on study-issued smartphones — the research-grade predecessor to today's CoughMonitor Suite

FEASIBILITY

Participants recorded a median of 23.9 hours per day of monitoring, showing that continuous, real-world cough capture was achievable across decentralized low-resource settings

The TB field is often slow to benefit from technological advances, but mobile health tools, like continuous cough monitoring, are an exciting avenue to make TB care more responsive and less burdensome on people with TB.

— Study Investigator, UCSF

WHAT THIS MEANS FOR YOUR STUDY

For sponsors and CROs developing TB and respiratory therapies, this study is early, peer-reviewed evidence that objective cough frequency could serve as a fast, decentralized signal of treatment response, potentially moving within days where sputum culture takes weeks. Proven here on Hyfe's research application across five countries, that same capability is now delivered on the wrist-worn CoughMonitor Suite platform, and is being carried forward in Hyfe's ongoing work with UCSF on cough and activity biomarkers as predictors of post-TB lung health.